

# Densitometer



Dolphin-Scan Plus is a professional Scanner for 2-D Gel, SDS-PAGE, Western Blot, X-ray film serving for high resolution image catch. Both transparency and reflectance mode are applicable for whether the samples are dry or wet.

## Features


- Transparency adapter for different film sizes
- Deliver the accurate images of gel sample
- Unique high speed blue light scanning
- Make batch scanning easily
- 1-D software available

## Applications

- Education, university, academic and research institute
- Analytical chemistry, molecular laboratory
- Food industrial and inspection laboratory
- Drug discovery, Biotechnology company
- Proteomics, genomics core facility
- Enzyme reaction, immunology

## 1D Software

- Support JPG, JPEG, BMP, TIFF, PCX, TGA, GIF format
- Comprehensive microtiter plate analysis
- Traceable image enhancement function
- Automatic data analysis with one click
- Specialized X-axis bands comparison
- Smiling gel identification ability
- Individual lane profile defining
- Iso-electric point simulation
- Intelligent colony counting
- Spot density calculation
- 3D image profile display



MW/Mass Determination

